

About this Catalog:

The online version of this catalog was designed with functional navigation to move between pages. Click the page numbers in the product list or the arrows to either side of the page numbers in the document to move to the products in the catalog.

Product List

Meatus and Incontinence Clamps	4
Control-Tip® Electrodes for rigid Cystoscopes	5
Flexible Electrodes for Rigid Cystoscopes	6
Electrodes for Flexible Ureteroscopes and Flexible Cystoscopes	7
Cutting and Coagulating Electrodes	8
Active Machine Plugs / Bi-Terminal Handles and Electrodes	9
Cords for Resectoscopes	10
Campbell Suprapubic Trocar and Cannula	11
Stamey Dorsal Vein Apical Retractor	11
Greenwald Urethral Suture Guides	12
Urethral Catheter Stylets	13
Urethral Sounds	14
Stamey Needles	14
Bougies a Boule	15

GREENWALD


Meatus and Incontinence Clamps







Cat. No	Description
U410	Whitver Meatus Clamp
U412	Hyams Meatus Clamp
U413	Zipser Meatus Clamp
U414R	Cunningham Incontinence Clamp, Regular
U414 J&L	Cunningham Incontinence Clamp, Juvenile & Large
U415	Baumrucker Incontinence Clamp


Control-Tip® Electrodes for rigid Cystoscopes

Features: Adjustable curve at electrode tip and 360° rotational feature result in precise tip control.
Color coded shaft insulation.
Precision machined nickel/silver electrode tips. Compatible with cat. no. M25 cord located on Page 9



			
RET 107A	RET 107B	RET 107F	RET 107K
Cat. No	Shaft Size	Description	
RET107A-05 RET107A-06 RET107A-07	5Fr. 6Fr. 7Fr.	Conical Tip	
RET107B-05 RET107B-06 RET107B-07	5Fr. 6Fr. 7Fr.	Bugbee Tip	
RET107F-05 RET107F-06 RET107F-07	5Fr. 6Fr. 7Fr.	Angular Tip	
RET107K-05 RET107K-06 RET107K-07	5Fr. 6Fr. 7Fr.	Ball Tip	
Insulation Color Code: 5Fr. - Light Gray, 6Fr. - Green, 7Fr. - Orange			

Flexible Electrodes for Rigid Cystoscopes



Features: Precision machined nickel/silver electrode tips. Color molded Teflon® shaft insulation. Compatible with cat. No. M25 cord on page 10



• **U51** Adaptor for CET Cystoscopic Electrodes. Permits use of hand-switch electrosurgical handle, eliminating the need for foot-switch.

CET 107



CET 107A



CET 107B



CET 107C



CET 107K



Cat. No	Shaft Size	Description
CET107-04 CET107-05 CET107-06 CET107-07 CET107-08	4Fr. 5Fr. 6Fr. 7Fr. 8Fr.	Pointed Tip
CET107A-04 CET107A-05 CET107A-06 CET107A-07 CET107A-08	4Fr. 5Fr. 6Fr. 7Fr. 8Fr.	Conical Tip
CET107B-03 CET107B-04 CET107B-05 CET107B-06 CET107B-07 CET107B-08	3Fr. 4Fr. 5Fr. 6Fr. 7Fr. 8Fr.	Bugbee Tip
CET107C-03 CET107C-04 CET107C-05 CET107C-06 CET107C-07 CET107C-08	3Fr. 4Fr. 5Fr. 6Fr. 7Fr. 8Fr.	Domed Tip
CET107K-04 CET107K-05 CET107K-06 CET107K-07 CET107K-08	4Fr. 5Fr. 6Fr. 7Fr. 8Fr.	Ball Tip

Shaft Insulation Color Code: 3Fr. - Black, 4Fr. - Red, 5Fr. - Gray, 6Fr. - Green, 7Fr. - Orange, 8Fr. - Blue

Electrodes for Flexible Ureteroscopes and Flexible Cystoscopes

Features: Precision machined nickel/silver electrode tips. Color coded Teflon shaft insulation. Compatible with cat no. M25 cord shown on page 10. Standard shaft lengths: 60, 90 and 120cm (custom lengths available on special order)



US1 Adaptor for EET Electrodes. Permits use of hand-switch electrocautery handle, eliminating the need for foot-switch.

*Teflon shaft insulation color code:
2Fr. - Blue • 3Fr. - Black • 4Fr. - Red • 5Fr. - Gray



Bugbee



Domed



Long Angular



Ball

Cat. No	Shaft Size	Description
EET107B-02 EET107B-03 EET107B-04 EET107B-05	2Fr. 3Fr. 4Fr. 5Fr.	Bugbee Tip
EET107C-02 EET107C-03 EET107C-04 EET107C-05	2Fr. 3Fr. 4Fr. 5Fr.	Domed Tip
EET107FC-03 EET107FC-04 EET107FC-05	3Fr. 4Fr. 5Fr.	Long Angular Tip
EET107K-04 EET107K-05	4Fr. 5Fr.	Ball Tip

Specify shaft length: (Unless specified otherwise, 120cm supplied on all 2 & 3Fr., 90cm supplied on all 4Fr. & 60cm on all 5Fr.)

Cutting and Coagulating Electrodes

Features: Nylon insulation. Precision machined 3/32" stainless steel shanks and tips (except SE79C which has a tungsten wire tip).
Fit most popular electrosurgical handles.

		
Cat. No	Description	
SE63	Turbinate Electrode	
SE79A	Fine Needle, 1"	
SE79C	Tungsten Round Loop, 1/4", 5/16", 3/8", 1/2", 3/4"	

Bi-Terminal Handle and Electrodes



Cat. No	Description
BT81	Bi-Terminal Handle with attached cord and two small banana machine plugs.
BT83	Hurd Turbinate Electrode

Active Machine Plugs

<p>M</p>	<p>Bovie Active Plug - Fits Aspen, Bovie AG, AS, CSV and 400 Series, Bard, E.M.S., Davol, Birtcher, Concept, 3M, Neomed, Valleylab</p>
<p>MH</p>	<p>Fits Birtcher Hyfreicator, National (Statham) Electricator</p>
<p>MSB</p>	<p>Fits Bovie Specialist</p>

Cords for Resectoscopes & Active Machine Plugs

Features: Machined Delrin instrument bushings. Lightweight, highly flexible, Silicon insulated conducting cable.
 *M25 cord is also compatible with electrodes on pages 5, 6 and 7.

M22



M22



M22



M22



M25



M25



M25



M25



Cat. No	Description (Standard Cord Length - 10ft)
M22	For Karl Storz #27050 series adult resectoscopes
M25	For Circon ACMI pre-USA adult resectoscopes **BEST SELLING CORD**

Note: Specify electrosurgical power unit used by make and model No., or select proper machine plug from those page 9.

Campbell Suprapubic Trocar and Cannula



Cat. No	Description
U440M	Campbell Suprapubic Trocar and Cannula

Stamey Dorsal Vein Apical Retractor

Features: Stainless steel throughout. Handle flange provides secure control and cephalad retraction. Retractor blades provide a wide space on both sides of the dorsal vein complex for placing a controlling suture.



Cat. No	Description
U530	Stamey Dorsal Vein Apical Retractor

Greenwald Urethral Suture Guides

Features: Precision machined, mirror-finish stainless steel components. Four suture guide grooves arranged at 90 degree intervals.



Cat. No	Description	Fr. Size
U501-24	Greenwald Urethral Suture Guide	24Fr.
U501-28	Greenwald Urethral Suture Guide	28Fr.

****Assist in reanastomosing urethra and bladder neck after radical prostatectomy****

Features: Three retractable gripping rods with non-slip ends. Three sets of paired suture guide grooves arranged at 120 degree intervals. Precision machined, mirror-finish stainless steel components.



Cat. No	Description	Fr. Size
U515-24	Greenwald No-Slip Tip Urethral Suture Guide	24Fr.
U515-28	Greenwald No-Slip Tip Urethral Suture Guide	28Fr.

Urethral Catheter Stylets

Features: Stainless steel throughout. Polished Finish.



Cat. No	Description
U320	Councill Stylet, 6Fr.
U321-04	Van Buren Curve Stylet, 4Fr.
U321-06	Van Buren Curve, Stylet, 6Fr.
U322-04	Malleable Tip Stylet, 4Fr.
U322-06	Malleable Tip Stylet, 6Fr.
U323	Guyon Curve Stylet, 6Fr.

Urethral Sounds



****Specify Fr. Sizes when ordering****

Cat. No	Description
U338	Hodgson Sound, Straight 8-14 Fr, 5" working length (Pediatric), 16-22 Fr, 7" working length (Adult). Short taper for periodic dilatation after urethrotomy
U340	McCrea Female Sound, 8-30Fr, 6" working length, 3/4" nominal bend height.
U343	Van Buren Sound, 10-30Fr, 9-1/2" working length, 1-3/4" nominal bend height.
U344	Walther Sound, 14-30Fr, 9-1/2" working length, 1-3/4" nominal bend height. Features proximal taper to relax the penile meatus

Stamey Needles



6 3/8" Working Length

Cat. No	Description
U403-0	Stamey Needle, Straight Tip
U403-15	Stamey Needle, 15° angled tip
U403-30	Stamey Needle, 30° angled tip

Bougies a Boule

Features: Mirror finish. Stainless Steel. Precision Machined.



****Specify Fr. Sizes when ordering****

Cat. No	Description
U419	Otis Bougie a Boule, Sizes 8-30 are stock sizes.
U420	Friedman-Otis Bougie a Boule, with 3 (graduated) olive tips on one stem. 12-14-16Fr., 18-20-22Fr., 22-24-26Fr.
U422	Thackston Bougie a Boule, with 2 (graduated) olive tips on each end. 14-16/18-20Fr., 22-24/26-28Fr.